



Packing instructions for gel electrophoresis box

The gel box contains fragile wires that can be damaged in transit if shipped with loose parts (gel trays, slats, and/or combs) inside the plastic box. Please follow these gel box packing guidelines when returning your MDLL kit.



Package the 2 gel slats and 6 gel trays together in bag.

Place bag of combs on side of gel box.



Place bag of gel tray ends on top of closed gel box.

Place bag gel trays and slats on side of gel box.



Gel box wires are platinum! Please do not pack anything inside the gel box where the wires are located.